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©Copyright owned by Qingdao ChuangMeng Instrument Co., Ltd

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Carefully read this User Manual to learn how to install and use the product correctly. After reading, properly keep the User Manual as a reference for future maintenance and repair.

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Address: No. 7 Wenzhou Road, City Northern District, Qingdao City, China

Production base: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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Manufactured by Qingdao ChuangMeng Instrument Technology Co., Ltd., 1103-model viscometer is the viscosity-measuring instrument for which the standard power adapter gives power supply, and it is applicable to 100~240V wide voltage belt, 50Hz or 60Hz power driving, and there are two speed types, namely 6-speed and 12-speed. See Diagram (I) for its structure and Table (I) for the details.

This is a typical coaxial cylinder viscometer in the annular space (shearing clearance) formed by placing the test fluid between the outer sleeve and the floater.

In viscosity measurement, when the outer sleeve rotates at a set speed, the fluid will produce a viscosity resistance, and thus

the viscosity of the fluid will be measured. When the outer sleeve produce a torque on the floater through the fluid, this torque acts on a precision spring, its deflection can be measured, and then make a comparison between the measuring result, measuring conditions and instrument constants. This system can imitate the real process in the industrial process.

Given a shearing rate, the coaxial cylinder rotational viscometer can make the viscosity measurement. This is a linear function relationship. For instance, the relational graph between the shearing stress and shearing rate is a straight line. In many cases, many fluids may not follow Newton's law, but their rheology is relatively close to Newton's law, so the viscometer can be still in use, and the viscosity can also be accurately measured. It is worth noticing that the calibration of the recommended 1103-model viscometer is a linear calibration following Newton's mode. This means that if the liquid sample does not conform to Newton's linear calculation method, and at this time, the viscometer will not be applicable. In this case, viscosity measurement and speed calculation will be subject to another linear calculation method, which is more applicable to the characteristics of such kinds of fluids.

When the standard outer sleeve, floater and torsional spring rotate in 300 rounds/minute, the measurement unit of the viscometer is centipoise (or milli pascal/second). At other rotating speeds, the read data must be conducted a series of simple multiplying conversion. The sixth part will present the calculation method of viscosity of plastic fluids (such as drilling fluid).

If choosing different rotating speeds or using different combinations of torsional spring-floater, the scope of shearing rate may have some changes. In order to expand the scope of shearing stress and thus measure all sorts of

liquids, we have designed different series of torsional springs, which can carry out reciprocal conversion easily.

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Laboratory technicians must be familiar with the operating procedures of the instrument and understand the potentially dangerous instruments and equipment. This instrument can be powered from 100 volts to 240 volts. Keep hands, clothes, and other objects away from the rotating parts of the instrument.

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The standard B1 floater equipped with 1103-model viscometer cannot be used to test the sample higher than 93°C.

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
There is a scale line at the 350ml of the slurry cup. Add the newly stirred test fluid to the liquid level line. The liquid level line on the outer sleeve can clearly show the proper immersion depth. It can be seen from

Diagram (I) that if the immersion depth exceeds the scale line, the bearing of the float may be damaged. If other containers are used, the distance between the bottom of the outer sleeve and the bottom of the container should be no less than 1.27 cm.

NOTE: If the standard BI floater is hollow, it cannot be used to test the sample higher than 93°C.

1103 model viscometer has six different speeds for testing, and the scope of speed is 3 rounds/minute to 600 rounds/minutes, and the rotating speed is decided by the control panel. When using it, first switch on the power, and then select the required speed button, read different values of shearing stress from the dial.

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Diagram (II) Control panel of 1103 model viscometer

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
The combination of R1-B1-F1 outer sleeve-floater-torsional spring is applicable to all viscometers. In order to calculate the shearing rate of the tested liquid, other types of combination of outer sleeve-floater-torsional spring may also be adopted. Some combinations may cause greater errors in the shearing stress readings, and such kinds of combinations do not meet the requirements.

##### A. 外套筒的拆除与安装 Assembling and disassembling of the outer cylinder

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Rotate the outer sleeve anticlockwise, and at this moment, the outer sleeve can be slowly removed, see Diagram (III). Rotate the outer sleeve clockwise, and make it rotate upward to the maximum position.


图（3）外套筒的拆装 Diagram (III) Assembling and disassembling of outer sleeve

B. 浮子的拆除与安装 Assembling and disassembling of the floater

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
The end side of the float shaft is cone-shaped, and it is inserted into a taper hole. When disassembling the floater, it requires to pull the floater downward and keeping it rotating anticlockwise at the same time. When assembling the floater, it requires to push the floater upward and keeping it rotating clockwise.

图（4）扭力弹簧的拆装 Table (VI) Assembling and disassembling of torsional spring

图（5）静负载校准装置  
Diagram (V) Static load calibration device

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Diagram (VI) Assembly drawing of torsion test

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1103 series viscometer should be inspected on a regular basis, and once a problem is found, it is necessary to carry out calibration and maintenance. Only when the instrument is calibrated regularly can the measurement precision be more accurate. Here two calibration methods are described to carry out calibration by applying torque to the float shaft.

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It is easily to operate static load calibration, and if the spring needs to be adjusted, its adjustment results can easily be confirmed. The standard liquid calibration requires a thorough calibration of the entire instrument. It is more accurate than static load calibration in checking the deflection of the floater, the eccentricity of the outer sleeve, etc. Please refer to Diagram (IV) B.

Carry out static load calibration with 1105 model torsional spring force-measuring and calibrating device, Please refer to Diagram (V).

1. Disassemble the outer sleeve, rotate the floater in the anti-clockwise direction and push it upward (the floater cooperates with the shaft taper), and then assemble the floater.
2. Take a piece of thin silk thread without elasticity, paste one end of the silk thread to the surface of the floater with a small piece of gummed paper, then rotate the silk thread towards the left for 2 ~ 3 rounds around the surface of the floater, and make the silk thread tie to the hook through a horizontal fast pulley or special force-measuring shelf.
3. Choose a weight based on (Table II).

砝码 Weight (g)	允许转动范围 Allowable rotation range (格 style)	实测 Measurement (格 style)	砝码 Weight (g)	允许转动范围 Allowable rotation range (格 style)	实测 Measurement (格 style)
5	21.55~22.21		40	172.45~177.70	
10	43.11~44.42		45	194.00~199.90	
15	64.66~66.63		50	215.57~222.13	
20	86.23~88.85		55	237.12~244.30	
25	107.78~111.06		60	258.68~266.56	
30	129.34~133.28		65	280.23~288.77	
35	150.89~155.49		70	301.80~310.99	

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This procedure is only applicable to the calibration of fluids that are in conformity with Newton's certification. The calibrated fluid can be 20,50,100,200 and 500cP. Every bottle of liquid that meets the ASTM standard is equipped with a comparison table of viscosity and temperature.

1. Ensure that the instrument to be tested is clean before immersing the outer sleeve and float into the standard liquid. If necessary, remove the outer sleeve and clean the float thoroughly. Make sure that the float shaft and outer sleeve are intact.

Attention: The batch number on the label of the standard liquid must match to the number on the viscosity / temperature chart.

2. Add the calibrated liquid to the liquid level line of the slurry cup, and place the slurry cup on the tray of the instrument. Lift the tray up until it is immersed to the proper depth. Refer to (Table I).
3. Put the thermometer in the sample to be tested, and choose a safe place to prevent it from breaking.
4. Start up, set 300 rounds for 3 minutes, and balance the temperatures of floater, outer sleeve, sample and the environment.
5. Record the readings on the dial for 300 rounds and 600 rounds, and the reading of the thermometer is exact to 0.1°C (0.15°F).

#### C. 扭力弹簧校准 Calibration of the torsional spring

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Please refer to Table (VI) for parts Identification.

Note: Make sure that the float shaft is not bent, and then start adjusting the torsional spring.

1. Disassemble the outer sleeve, rotate the floater in the anti-clockwise direction and push it upward (the floater cooperates with the shaft taper), and then assemble the floater.
2. Take a piece of thin silk thread without elasticity, paste one end of the silk thread to the surface of the floater with a small piece of gummed paper, then rotate the silk thread towards the left for 2 ~ 3 rounds around the surface of the floater, and make the silk thread tie to the hook through a horizontal fast pulley or special force-measuring shelf.
3. Hang a 5~65g standard weight for calibration, and the dial readings should conform to Testing Table of Stiffness and Linearity of Torsional Spring

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When operating at 300rpm, the reading on the R1 B1 F1 combined dial is Newtonian viscosity. If using other springs, it is necessary to multiply the dial readings with “f” factor (spring constant).

Determine Newtonian viscosity cP with the viscometer, and apply the following formula:

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Example: use a R2 B1 combination at the speed of 600 rpm and f5.0 spring, one dial reading is 189.

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Note: Calibration-intended standard liquid has the scope of viscosity and the standard R-B-F combination will be used for testing. If choosing improperly, the test data will be affected.

外套筒-浮子组合 Outer sleeve- floater combination	外套筒-浮子因子 Outer sleeve- floater factor
R1-B1	1.000
R1-B2	8.915
R1-B3	25.392
R1-B4	50.787
R2-B1	0.315
R2-B2	8.229
R2-B3	24.707
R2-B4	49.412
R3-B1	4.517
R3-B2	12.431
R3-B3	28.909
R3-B4	57.815

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Example: the required setting is F1 torsional spring, and the constant of torsional spring is the combination of 386 Dynes/cm/° def and R1 B1. Use the 50g weight, and the formula is:

$$K1 = 50 \times 1.725 \times 981 / 386 = 219.2$$

## 3. 数据测试及计算 Data testing and calculation:

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Adjust the room temperature at 20±5°C, and work in strict accordance with the "test operation method". If making measurement at the well site, reduce the time delayed for taking the sample as much as possible, and record the sampling site and condition on the measurement table.

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Adjust the speed of the instrument to 300r/min, until the readings on the dial are constant, and the readings are absolute viscosity values.

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1) Adjust the speed of the instrument to 600r/min, until the readings on the dial are constant, and 1/2 of the readings are apparent viscosity values.

2) Adjust the speed of the instrument to 300r/min, and the differences between its readings and the readings of 600r/min are apparent viscosity values.

3) Stir the drilling fluid at a high speed for 10 seconds, and after rotating at the speed of 3r/min, the maximum value is the initial shearing force. Set it aside for 10 minutes, and record static shearing force.

$$\text{视粘度: } \eta_{\text{视}} = 1/2 \times 600 r/min (\text{读数}) \quad \text{mPa} \cdot \text{s}$$

$$\text{Apparent viscosity: } \eta = 1/2 \times 600 r/min (\text{reading}) \quad \text{mPa} \cdot \text{s}$$

$$\text{塑性粘度: } \eta_{\text{塑}} = 600 r/min (\text{读数}) - 300r/min (\text{读数}) \quad \text{mPa} \cdot \text{s}$$

$$\text{Plastic viscosity: } \eta_{\text{plastic}} = 600 r/min (\text{reading}) - 300r/min (\text{reading}) \quad \text{mPa} \cdot \text{s}$$

$$\text{动切力: } \tau_0 = 5.11 (300r/min \text{ 读数} - \eta_{\text{塑}}) \quad \text{Pa}$$

Dynamic shearing force:  $\tau_0 = 5.11$  (300r/min reading  $-\eta$  plastic) Pa

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$\tau$  (initial) =  $5.11 \times 3r/\text{min}$  (reading) Pa (Set it aside for 1 minute)

$\tau_{终} = 5.11 \times 3r/\text{min}$  (读数) Pa (静置 10 分钟)

$\tau$  (final) =  $5.11 \times 3r/\text{min}$  (reading) Pa (Set it aside for 10 minutes)

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Its flowing characteristic is that it begins to flow when there is shearing stress, but the viscosity decreases with the increase of the shearing stress, and the flow of the pseudoplastic fluid obeys the power function. The formula is:

$$\tau = k \left( \frac{dv}{dx} \right)^n \quad \lg \tau = \lg k + n \lg \frac{dv}{dx}$$

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n——epidemic index its value is between 0 and 1

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## VIII.故障排除与维护 Trouble-shooting and maintenance

### 1.故障排除 Trouble-shooting

故障 Fault	原因 Cause
刻度盘读数不稳 The dial readings are unstable	1.浮子轴轴承生锈 1.The bearing of the float shaft is rusted 2.浮子轴弯曲 2.The float shaft is bent 3.外套筒失准 3. The outer sleeve is out of alignment
数据不准 The data is inaccurate	1.浮子轴轴承生锈 1.The bearing of the float shaft is rusted 2.浮子轴弯曲 2.The float shaft is bent 3.外套筒弯曲 3. The outer sleeve is bent 4.扭力弹簧损坏或安装不正确 4.The torsional spring is damaged or its installation is incorrect 5.电机需要更换 5. The motor needs to be replaced
噪音过大 The noise is much too loud	1.电机故障 1. Motor failure 2.壳体螺丝松动或安装不当 2. The shell screw is loosen or the installation is improper.
外套筒径跳过大 The radial beat of the outer sleeve is too large	1.外套筒损坏 1. The outer sleeve is damaged 2.传动轮或传动皮带损坏 2. The driving wheel or driving belt is damaged.
按键失灵 Keypad failure	1.主控制板故障 1. The main control panel is damaged
电机不运行 The motor is not operated	1.电机损坏 1. The motor is damaged 2.驱动器损坏 2. The driver is damaged 3.电源插头未插好 3. The attaching plug is not well inserted

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- Promptly clean the floater and the outer sleeve after each test, and regularly inspect the indentations, wear and tear or other damages.
- The normally used viscometer does not require oiling or lubrication.

c. Take down the floater and the outer sleeve in the transporting process to avoid that the float shaft is bent and the outer sleeve is damaged, and test the bearing of the float shaft on a regular basis.

d. In the absence of a sample, if the operating instrument observes the operation of the floater and outer sleeve at 3 rpm or 6 rpm, there should not be the fluctuations of over  $\pm 1$ . The insensitive float shaft should be replaced in a timely manner. The dials should be maintained by professional maintenance personnel.

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Production address: No. 3 Xinghai Road, Liuting Street, Chengyang District, Qingdao

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## Qingdao Chuangmeng Instrument Co., Ltd. Packing list
Manufacturing enterprise: Qingdao Chuangmeng Instrument Co.,Ltd.